

PACKING GUIDE FOR SHIPPING PACKAGES

Correct packaging is the best guarantee for the package to reach its destination intact and with the correct temperatures.

LabAnalysis offers Customers everything they need to prepare the package for the sample shipment.

Operations to be performed are the following:

CONTAINER INTEGRITY

Always check container integrity. The container must be intact, sturdy, without holes, without scratches or damages of any kind. The goods will be more protected in a resistant box and when you receive it, you can visually check the treatment of your package during the journey!

INTERNAL PROTECTION

Create an inner layer with some insulating material such as polystyrene or bubble wrap, making sure you don't leave any gaps, so that the goods inside don't move. Insert a plentiful quantity of synthetic ice (icing box, cooling plates).

FIRM CLOSURE AND CLEAN SURFACES

The package must be securely sealed with packing tape. Make sure that the package is difficult to open by applying adhesive tape on all ends of the package and on all junctures.

Delete any labels that are already present on the package and that are not part of the shipment. Therefore, delete any reference not related to the shipment

